


# 100% Silicone Foley Catheter

## Standard, 2-Way


**Intended use:**

The 100% Silicone Foley Catheters are indicated for transurethral drainage of urine from the bladder. The adult 2-way catheters can be retained in the human body for up to 90 days. The standard 2-way catheters (sizes 12 to 18) can also be used for replacement of suprapubic catheters for drainage of the bladder (exchange only). The products are intended for single use only.


**Product range:**
**Silicone Foley Catheter, Standard, 2-Way, 10ml balloon**

	REF*	Total Length, mm (±10)	Balloon Capacity, ml	CH/FG	Tip	For Suprapubic Exchange
 UM41151410	UM41151210	400	10	12	Nelaton	Yes
	UM41151410	400	10	14	Nelaton	Yes
	UM41151610	400	10	16	Nelaton	Yes
	UM41151810	400	10	18	Nelaton	Yes
	UM41152010	400	10	20	Nelaton	-
	UM41152210	400	10	22	Nelaton	-
	UM41152410	400	10	24	Nelaton	-
	UM41152610	400	10	26	Nelaton	-

**Silicone Foley Catheter, Standard, 2-Way, 30ml balloon**

 UM41151830	UM41151630	400	30	16	Nelaton	-
	UM41151830	400	30	18	Nelaton	-
	UM41152030	400	30	20	Nelaton	-
	UM41152230	400	30	22	Nelaton	-
	UM41152430	400	30	24	Nelaton	-
	UM41152630	400	30	26	Nelaton	-

\* Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

# 100% Silicone Foley Catheter

Standard, 2-Way



## Material specifications:

**Contains latex?** . . . . . No  
**Contains phthalates?** . . . . . No  
**Catheter shaft and funnel** . . . . . Silicone  
**X-ray line** . . . . . Silicone  
**Balloon** . . . . . Silicone  
**Tip** . . . . . Silicone, Barium sulphate  
**Hard valve** . . . . . Polyvinyl chloride, Polypropylene,  
 Acrylonitrile-butadiene-styrene, Silicone, Stainless steel  
**Sleeve** . . . . . Polypropylene

## Certifications:

**CE-marking?** . . . . . Yes (0123)  
**Product classification** . . . . . Class IIB (EU)  
**Quality system certified to** . . . . . ISO 13485  
**Global Medical Device Nomenclature (GMDN)** . . . . . 34917  
 For suprapubic exchange (sizes 12 to 18): 34924

## Sterilization:

**Sterile?** . . . . . Yes  
**Sterilized by** . . . . . Ethylene oxide  
**Shelf life** . . . . . 5 years

## Packaging specifications:

REF	Primary package/ Single peelpack		Secondary package/ Shelf box		Tertiary package/ Shipper				Pallet
	Qty, pcs	Dimensions, mm (LxW)	Qty, pcs	Dimensions, mm (LxWxD)	Qty, pcs	Dimensions, mm (LxWxD)	Gross Wt. kg	Volume m <sup>3</sup>	Qty, pcs
UM41151210	1	500 x 66	10	510 x 80 x 45	200	520 x 340 x 250	6.0	0.044	4,800
UM41151410	1	500 x 66	10	510 x 80 x 45	200	520 x 340 x 250	6.0	0.044	4,800
UM41151610	1	500 x 66	10	510 x 80 x 45	200	520 x 340 x 250	6.0	0.044	4,800
UM41151810	1	500 x 66	10	510 x 80 x 45	200	520 x 340 x 250	6.0	0.044	4,800
UM41152010	1	500 x 66	10	510 x 80 x 45	200	520 x 340 x 250	6.0	0.044	4,800
UM41152210	1	500 x 66	10	510 x 80 x 45	200	520 x 340 x 250	6.0	0.044	4,800
UM41152410	1	500 x 66	10	510 x 80 x 45	200	520 x 340 x 250	6.0	0.044	4,800
UM41152610	1	500 x 66	10	510 x 80 x 45	200	520 x 340 x 250	6.0	0.044	4,800
UM41151630	1	500 x 66	5	510 x 80 x 35	250	520 x 430 x 375	7.0	0.084	1,500
UM41151830	1	500 x 66	5	510 x 80 x 35	250	520 x 430 x 375	7.0	0.084	1,500
UM41152030	1	500 x 66	5	510 x 80 x 35	250	520 x 430 x 375	7.0	0.084	1,500
UM41152230	1	500 x 66	5	510 x 80 x 35	250	520 x 430 x 375	7.0	0.084	1,500
UM41152430	1	500 x 66	5	510 x 80 x 35	250	520 x 430 x 375	7.0	0.084	1,500
UM41152630	1	500 x 66	5	510 x 80 x 35	250	520 x 430 x 375	7.0	0.084	1,500

## EAN codes:

REF	Primary package	Secondary package	Tertiary package
UM41151210	(01)30768455149088	(01)50768455149082	(01)70768455149086
UM41151410	(01)30768455149095	(01)50768455149099	(01)70768455149093
UM41151610	(01)30768455149101	(01)50768455149105	(01)70768455149109
UM41151810	(01)30768455149118	(01)50768455149112	(01)70768455149116
UM41152010	(01)30768455149125	(01)50768455149129	(01)70768455149123
UM41152210	(01)30768455149132	(01)50768455149136	(01)70768455149130
UM41152410	(01)30768455149149	(01)50768455149143	(01)70768455149147
UM41152610	(01)30768455149156	(01)50768455149150	(01)70768455149154
UM41151630	(01)30768455149163	(01)50768455149167	(01)70768455149161
UM41151830	(01)30768455149170	(01)50768455149174	(01)70768455149178
UM41152030	(01)30768455149187	(01)50768455149181	(01)70768455149185
UM41152230	(01)30768455149194	(01)50768455149198	(01)70768455149192
UM41152430	(01)30768455149200	(01)50768455149204	(01)70768455149208
UM41152630	(01)30768455149217	(01)50768455149211	(01)70768455149215



# 100% Silicone Foley Catheter





Standard, 2-Way



## Packaging material:

Primary package	Sleeve - Ethylene-vinyl acetate, Peelpack pouch - Coated paper, Film
Secondary package	Ivory cardboard
Tertiary package	Corrugated board

## Storage and handling:

	<b>FRAGILE</b> - Contents of the transport package are fragile and shall therefore be handled with care
	<b>KEEP DRY</b>
	<b>THIS WAY UP</b> - Indicates correct upright position of the transport package
	<b>DO NOT USE HAND HOOKS</b> - Hooks are prohibited for handling the transport package



**Well Lead Medical Co., Ltd.**

C-4 Jinhu Industrial Estate  
Hualong, Panyu, Guangzhou  
511434 P.R. China

Made in China



**Shanghai Intl. Holding Corp. GmbH**

Eiffestrasse 80, 20537 Hamburg  
Germany

Distributed by



ConvaTec Limited  
First Avenue, Deeside Industrial Park  
Deeside, Flintshire, CH5 2NU, UK

[www.convatec.com](http://www.convatec.com)

**For further information on these products, please contact your local ConvaTec sales representative.**

With reservations for misprints and alterations in product specifications. Printed version not controlled.

All trade marks and logos are the property of a ConvaTec group company.

© 2018 ConvaTec Inc. AP-019414-MM, ver. OCT 2018